English
Reprocessing instructions for reusable WAL cannulae

This document provides detailed instructions on cleaning, disinfection and sterilization of the WAL infiltration and irrigation/ aspiration cannulae (including *harvest* cannulae, Biofill cannulae) produced by Human Med AG.

Cannulae that are supplied **non-sterile must** be thoroughly cleaned and sterilized before use.

Warning notices:	WAL cannulae must be reprocessed exclusively in automated procedures! The use of ultrasonic baths for processing cannulae with anodized aluminum handles may affect the quality of the anodized surfaces if the cannulae touch each other. Devices with anodized surfaces are damaged by alkaline (pH > 10) cleaning agents and disinfectants. Therefore, only cleaning agents recommended for anodized aluminum should be used. Please consult the manufacturer of the instrument cleaning agent or disinfectant in this regard. Cleaning, disinfection and sterilization shall only be carried out according to the state of the art by trained staff in facilities specifically intended for this purpose. Note: When processing medical devices take care in working carefully according to the local safety regulations.
Information on clinical use:	When using alkaline cleaners for processing cannulae with anodized aluminum handles, it is to be expected that the color of the cannula handles changes and the anodized coating peels off. If cannulae with anodized aluminum handles shall be reprocessed with alkaline cleaners and afterwards used for liposuctions with subsequent autologous fat transfer or fat collection for laboratory purposes, no unlimited guarantee is given for the biocompatibility of the collected fat. Only cannulae with stainless steel handles are to be used for such surgical interventions. These cannulae were successfully tested by Human Med with regard to microbial safety, and biocompatibility; they can safely be reprocessed with alkaline cleaners.
Limitations on reprocessing:	Provided that the products are used as intended by the manufacturer, the end of their service life is determined by normal wear and tear and use of the appropriate cleaning agents. Frequent reprocessing may affect the outside appearance of the WAL cannulae. Prior to each reprocessing cycle the cannulae must be checked for damages. Do not use damaged cannulae.
REPROCESSING PROCEDURE	
Instructions:	
Preparation and transport::	Remove surface contamination using a disposable cloth or paper towel. Take care that blood and tissue residues do not dry on the products. We recommend reprocessing immediately after use. Do not use fixing agents or hot water (> 40°C, 104°F) since these will result in the coagulation and fixation of proteins to the instruments' surfaces, affecting the success of subsequent cleaning procedures. Special care must be taken during transport to ensure that the distal ends of the WAL cannulae are not damaged. They should be transported in a closed container to avoid contamination to the environment. You are recommended to start the reprocessing of the WAL cannulae immediately after use.
General information:	The WAL cannulae with anodized aluminum handles must exclusively be treated with enzymatic cleaning agents. WAL cannulae with stainless steel handles may also be treated with alkaline cleaners. We recommend using an automated washer-disinfector for MIS instruments and the processing program MIS INSTR. NEUTRAL. Note: When processing medical devices always handle with care in accordance with local Health & Safety procedures.
Pre-cleaning:	Pre-clean the interior of the cannula handle using a pipette brush (Ø ca. 20 mm), then shortly flush the cannula handle by hand in a bowl filled with deionized water. Dispose of or reprocess the pipette brush after use.
The following steps must be performed in a validated process using a qualified washer-disinfector:	
Use of a rinse adapter:	Insert the rinse adapter (we recommend special rinsing lances produced by 3 Mach GmbH) into the cannula handle and connect the rinsing tube of the washer-disinfector to the rinse adapter by attaching the Luer adapter of the rinsing tube to the Luer connector of the rinse adapter. The rinse adapter is to be reprocessed according to the same procedure as the WAL cannulae.
Cleaning:	Pre-rinse with cold water in the washer-disinfector for 15 minutes. The cleaning cycle is done in 55°C hot water, mixed with an enzymatic cleaning agent (neodisher MediZym or Thermostept RKN-zym are recommended; concentration as indicated by the manufacturer of the cleaning agent) for a temperature holding time of at least ten (10) minutes. The use of alkaline cleaning agents is also approved for the cleaning of WAL cannulae with stainless steel handles. Afterwards flush with a 0,1 % concentrated acid agent (Neodisher Z or Thermostept NKZ are recommended) diluted in 20°C warm water for ten (10) minutes.
Disinfection:	Rinse with 20°C warm, deionized water for at least five (5) minutes. Thermal disinfection is done with softened water. Human Med carried out the validated reprocessing at a temperature of 93°C and at a temperature holding time of at least ten (10) minutes. You are recommended to validate the temperatures for disinfection on your own according to the AO concept (see standard EN DIN ISO 15883).
Drying:	Dry at a temperature of 60°C for at least 15 minutes.
Cleaning check:	Check all surfaces, holes and lumens for visible contamination. Contaminated products must be returned to the reprocessing cycle.
Maintenance:	Apply a small quantity of surgical grade silicone spray lubricant to the area of the handle where the WAL cannula is attached.
Visual check:	Visually check the cannulae for damage and wear (e. g. firm seating of the suction tube in the cannula handle and intactness of the cannula's distal end). Reject any worn or damaged cannula. Note: If any cannula has to be returned to Human Med or your supplier, it must be cleaned, disinfected and sterilized before. Enclose the relevant certificates stating this.
Packaging:	In accordance with standard ISO 11607 the cannulae must be individually packed in sterile packaging. The packaging must be large enough so that the sealing is neither subject to stress nor can be damaged.
Sterilization:	Steam sterilization using a fractionated vacuum process at a temperature of 121°C (250 F) for at least 15 minutes, and steam sterilization at 134°C (273°F) for at least three (3) minutes have been validated. Prior to each sterilization process check the batch process records for the respective sterilizer to ensure that the required parameters are being complied with.
Storage:	Store in a closed cabinet, protected from dust and humidity and extreme fluctuations in temperature. The shelf life is determined by the product specifications of the sterile packaging used.

The instructions outlined above have been validated by Human Med as being appropriate for the preparation of the WAL cannulae for reuse.

Information about validation of reprocessing

The following test instructions, materials and machines have been used for validation:

Cleaning agent (automated process): Only for cannulae with stainless steel handles: Neutralizer (automated process): Washer-disinfector:	Neodisher MediZym, Dr. Weigert mediclean forte, Dr. Weigert Neodisher Z, Dr. Weigert GETINGE, CM 305, program MIS INSTR. NEUTRAL
--	---

Contact Human Med AG for detailed information on the validations.

If the above described chemicals and machines are not available, it is the user's responsibility to validate the reprocessing method used.

It is the responsibility of the reprocessor to ensure that the reprocessing actually performed with the used equipment, materials and personnel in the reprocessing facility achieves the desired results.

This requires validation and routine monitoring of the process. Likewise any deviation from the instructions provided must be properly evaluated by the reprocessor for its effectiveness and potential adverse consequences.

For further information please contact Human Med under +49 (0) 385 395700 or the supplier of the medical device.

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I dettagli concernenti la validazione possono essere richiesti alla ditta Human Med AG.
Vyše uvedené pokyny byly firmou Human Med AG validovány jako vyhovující pro WAL kanyly.
Podrobné údaje o validaci si můžete vyžádat od firmy Human Med AG.