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CE 2797
www.merit.com

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222

402533002MLP_001 ID 2020-07-29

EN

Intended use:

The devices contained in the face lipofilling kit are used to perform autologous free fat transfer. The procedure usually involves harvesting adipose tissue with liposuction, preparing the aspirated tissue with various techniques and injecting the prepared tissue into the receptor sites.

Target population:

For human use - Do not use on children or pregnant women

Intended users:

Health professionals only

Contraindications and side effects:

Contraindications and side effects of infiltration and sampling cannulas are described in the NUP-LT leaflet 'TAKING AND INJECTION CANNULAS FOR LIPOTRANSFER' attached to the kit

The contraindications and side effects of filling cannulas are described in the NUP-MI leaflet 'MICRO CANNULAS AND FILLER INJECTION KITS' attached to the kit

Instructions For Use:

The collection and injection of autologous adipose tissue are performed in 1 surgical step under general or local anesthesia.

An anesthetic solution is injected subcutaneously at the planned insertion points of the cannulas. This injection is performed using a needle and syringe not provided in the kit. The incision site is disinfected and a small incision is made using a 14G needle, ref IN-21X40 or 8300013707. Small incisions are made in different places depending on the area to be collected but are always designed to hide the surgical scar.

Tumescence of the area to be sampled can be performed according to the user's needs. The infiltration solution is infused into the tissues before liposuction using an infiltration cannula (not supplied in the kit).

The fat is collected using the 8-hole $\varnothing$ 2.1 x 150 mm fat harvesting cannula, ref CP141508H connected to a 30 ml VacLoK AT® syringe, ref VAC230. Adipose tissue can be prepared by centrifugation after being transferred to 10 ml syringes (not supplied in the kit) using the FLL / FLL connector, ref FD0993922-T. Once filled, the syringes are sealed with a Luer Lock cap, ref FD-0001711-MC before centrifugation.

The adipose tissue is then transferred into smaller syringes (not supplied in the kit) using the FLL / FLL connector, ref FD0993922-T in order to allow reinjection using a subcutaneous filling cannula.

The injection site is disinfected and small incisions are made using an 18G needle, ref IN-12X40 or 4790012040 which will allow the insertion of the 18G reinjection cannula.

The adipose tissue is thus injected via a cannula for subcutaneous filling $\varnothing$ 1.25 x 110 mm, ref MI18110KR or $\varnothing$ 0.81 x 40 mm, ref MI2140KR or $\varnothing$ 0.81 x 70 mm, ref MI2170KR depending on the area injected.

Warnings and precautions:

The warnings and precautions of the infiltration and harvesting cannulas are described in the NUP-LT leaflet 'HARVESTING AND INJECTION CANNULAS FOR LIPOTRANSFER' attached to the kit

The warnings and precautions for filling cannulas are described in the NUP-MI 'MICRO CANNULA AND FILLER INJECTION KITS' notice attached to the kit.

CAUTIONS OF THE DEVICES CONTAINED IN THE KIT:

Read the instructions and precautions on each device before use.

Cautions for the needle ZARYS dispoFINE, Reference IN-21X40 and IN-12X40:

Device designed for intravenous, intramuscular and hypodermic injections.

- Before using the product, inspect the packaging to make sure it is not damaged.
- Open the packaging and pull out the needle
- Attach the needle to a syringe with luer or luer-lock tip
- Carefully remove the protective cap
- Perform an injection according to the effective aseptic medical techniques
- Discard the needle into a sharps container after the injection
- If the needle is damaged or contaminated, replace it with a new one before you perform an injection

If this product is to be used during drug preparation or aspiration from a vial or similar container, it is recommended to use a needle with the diameter not exceeding 0,8 mm (no. 8)

Considering the risk of contaminating the drug with pieces of a vial stopper, it is recommended to perform drug aspiration and preparation with devices designed specifically for drug aspiration.

Precautions concerning the Vac Lok AT™ syringe, reference VAC230:

English

VAC ▶ LOK ▶ AT™

Vacuum Pressure Syringe

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION
The VacLok AT™ Vacuum Pressure Syringe, is a general piston syringe constructed using a barrel, plunger, piston seal and locking mechanism. The plunger may be locked in any position along the length of the barrel to maintain a constant vacuum pressure.

INDICATIONS FOR USE
The VacLok AT Vacuum Syringe is used to inject fluids into, or withdraw fluids from the body. It can also be used in cases where a vacuum syringe is preferred (e.g. thrombus, abscess fluid, bile, urine etc.).

CLINICAL BENEFITS
Administration or withdrawal of fluids or the creation of a vacuum during appropriate procedures.

PRECAUTIONS

- Read manufacturer's instructions prior to use
- Rx Only Caution: Federal (USA) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of the device.
- Inspect the device prior to use to verify that no damage has occurred during shipping.
- Do not use after expiration date

REUSE PRECAUTION STATEMENT
For single patient use only. Do not reuse, reprocess or re sterilize. Reuse, reprocessing or re sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re sterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient.

CONTRAINDICATIONS
None known

INSTRUCTIONS FOR USE

- Using aseptic technique, open the package containing the sterile product.
- Inspect for damage or improper assembly.
- To activate vacuum pressure (Figure 1):
 - Pull plunger back to desired location.
 - Twist plunger clockwise until it stops and locks into place.
- To disconnect vacuum pressure (Figure 2):
 - Twist plunger counter-clockwise until plunger stops.
- To aspirate fluids without locking the syringe, pull plunger back and maintain pressure without turning the plunger.


Figure 1


Figure 2

Nominal capacity accuracy per ISO 7886-1

	Single use
	Caution
	Non-pyrogenic
	Do Not Use if Package is Damaged and Consult Instruction for Use
	Manufacturer
	Use by date: YYYY-MM-DD
	Authorized Representative in European Community
	Catalog number
	Lot number
	Date of Manufacture: YYYY-MM-DD
	Medical Device
	Unique Device Identifier
	Consult Instructions for Use For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID. For printed copy, call U.S.A. or EU Customer Service
	Do not re sterilize
	Single sterile barrier system
	Sterilized Using Ethylene Oxide

CE 2797
www.merit.com

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222

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