



**MANUFACTURER:** TULIP MEDICAL PRODUCTS

**PRODUCT DESCRIPTION:** Anaerobic Transfer Luer Lock to Luer Lock. Supplied Sterile.

**MODEL NAME:** GEMS ANAEROBIC TRANSFER LUER LOCK TO LUER LOCK

**MODEL DESCRIPTION:** GEMS Anaerobic Transfer, Single-use, for Adipose Tissue.

**PRODUCT CODES AND SIZES:**

Single-Use Anaerobic Transfer, Luer Lock to Luer Lock	DATLLLL 1,2 MM
Single-Use Anaerobic Transfer, Luer Lock to Luer Lock	DATLLLL 1,4 MM
Single-Use Anaerobic Transfer, Luer Lock to Luer Lock	DATLLLL 2,4 MM

**PRINCIPAL MATERIAL COMPONENTS:**

- ABS Transparent

**SUPPLY CONDITIONS AND CONSERVATION:**

- Packaged Sterile, Single-Use only

**STERILIZATION AND EXPIRATION OF PRODUCT:**

Sterile:	Yes	Re-sterilizable:	No
Expiration:	5 years	Single-use:	Yes

**STERILIZATION METHOD:** Do Not Resterilize

**CE MARK:** SELF DECLARED CLASS 1 CE MARKED

**LATEX FREE:** YES

**CE Conformity Assessment:** Class IIA, Rule 6, Annex II of MDD 93/42/EEC Council Directive

**CND CODE:** A0701

**INSTRUCTIONS FOR USE AND WARNINGS:**

SEE PACKAGE INSERT